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

# SURFACE AREA/POROSITY ANALYSIS

## GAS SORPTION ANALYSIS OF SOLIDS

Gas sorption analysis using nitrogen, argon or krypton enables the determination of specific surface areas and pore size distributions. The ASAP2420 (by Micromeritics), available at the University of Salzburg, can be used for surface area and porosity analysis of porous solids, such as powders or monolithic samples.

N<sub>2</sub>-sorption isotherm of a micro-, meso- and macroporous silica material.

### APPLICATIONS:

Types of gas sorption analysis measurements typically performed by ASAP2420 are for example:

- micropore analysis
- surface area calculations using BET-method
- determination of pore size and pore shape
- pore size distributions of mesopores
- determination of particle size
- characterization of surface properties of metal powders, ceramic materials, catalysts, fillers, building materials, pigments and many other solids.

### METHOD:

ASAP2420 enables the automatic determination of physisorption isotherms of adsorbed gas molecules at porous surfaces under isothermal conditions (77 K for N<sub>2</sub>). Adsorption data, such as quantity adsorbed at certain relative pressures, can be used to evaluate the specific surface area, pore size and pore size distribution, pore shape and surface properties (hydrophobic or hydrophilic) of the samples by different models (BET, BJH, DFT, etc.).

### BENEFITS:

Small sample amount (mg), short measuring times, high throughput.

### SPECIFICATIONS:

Degas System: 12 degas ports  
Temperature range: 25-450°C  
Evacuation rate: 1-50 mmHg/s

### Analysis System:

6 analysis ports  
5 analysis ports for micropore measurements  
Isothermal jackets for sample- and reference tube

### SERVICE:

Samples can be measured within 4 to 5 days.

### CONTACT:

Technology Transfer  
Dr. Erika Hebenstreit LL.M.  
Forschungsservice  
Kapitelgasse 4-6  
5020 Salzburg

Tel.: +43-(0)662-8044 2451  
Erika.Hebenstreit@sbg.ac.at  
www.uni-salzburg.at/forschungsservice/techtransfer

### GAS SORPTION:

Univ.-Prof. Dr. Nicola Hüsing  
Department of Materials  
Science and Physics  
Hellbrunnerstraße 34  
5020 Salzburg

Tel.: +43 (0)662-8044 5404  
nicola.huesing@sbg.ac.at  
www.uni-salzburg.at/mw